

Paging Gateway

Providing centralized, effective and easy to manage paging services in today's organizations is a challenging task for any communications manager. Controlling cost, providing enhanced and efficient paging services while maintaining control over pager ID numbers, on site transmitters and Radio Common Carriers (RCC's) has become a full time job.

Amcom Software's Paging Gateway was developed with your organization in mind. Paging Gateway is a seamlessly integrated paging platform that enables staff to generate pages – voice, numeric or alphanumeric – with or without operator intervention.

Integrated Solution

- Intimately integrates with other Amcom products, thus eliminating the need for additional databases to maintain

Standardize Paging

- Whether sending alphanumeric pages to an RCC, or numeric and voice pages to your on-site transmitter, Paging Gateway makes this detail transparent and consistent with both the paging party and the pager carrier

Page Forwarding and Coverage

- Allows end-users to forward pages to other subscribers or other pagers under password control both on demand and on a scheduled basis

Fully Integrated Systems

- Whether you maintain an on-site transmitter or rely heavily on an RCC – or both – Paging Gateway can make this transparent to both the party being paged and the end-user carrying the pager. Use of separate pager ID's allows administrators to change out paging providers without retaining end users on paging access or new pager numbers

Traditional Operator Functions Which Eliminate Costs

- Paging Gateway accepts inputs from other Amcom modules and allows IVR sessions for numeric and voice pages, and accepts TAP or IXO input from other alpha paging systems or nurse call systems. Paging Gateway connects to PURC-compatible and other on site transmitters without your organization having to purchase additional paging encoders



For more information about Paging Gateway or Amcom Software, call 888-337-9274 or visit www.amcomsoftware.com.